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## Your Advantages

- up to 1.12 m<sup>2</sup> of workspace
- modular construction
- space-saving set up directly at the wall
- easy opening and closing of the sliding door
- corrosion-resistant, plastic laminated construction
- adjustable flow rate
- flow indication via LED
- connections for central building control systems
- multiple configuration

The laminar flow products from asecos enable you to easily establish a cleanroom workplace with particle reduced air in a normally ventilated area. Perfect for all areas where a particle free environment is required during processing, preparation or storage.

### The Laminar Flow Box Technology

The EC fan of the asecos Laminar Flow Box FBS sucks in the ambient air and extrudes it through high performance filters. This generates a laminar flow (air flows in parallel ways without any whirls) which captures the particles inside of the Laminar Flow Box and guides them through a perforated metal plate to the outside. Inside the Laminar Flow Box is now a cleanroom, type class 5 in accordance with EN ISO 14644-1.

### The Filters

The asecos Laminar Flow Box FBS uses a filter type H14 with laminatisator. This filter has a filtration efficiency of 99.995 %. By the use of this filter the Laminar Flow Box provides an isolation factor of 104. Air quality inside the box is improved by 10.000 times.

Diagram: Laminar Flow Box FBS Series

### The Product Range

The asecos Laminar Flow Box FBS product range offers solutions for varied tasks. Additional configurations complete the modular construction, e.g interior lighting, ceiling mount or base frames.

#### Laminar Flow Modules FMS-Series – filter module without box

##### Dimensions filter module without box - FMS-Series:

| Type                        | (W x D mm) | kg | Order No. |
|-----------------------------|------------|----|-----------|
| Laminar Flow Module FMS 24  | 610 x 400  | 20 | 21731     |
| Laminar Flow Module FMS 37  | 610 x 610  | 23 | 21733     |
| Laminar Flow Module FMS 56  | 915 x 610  | 30 | 21735     |
| Laminar Flow Module FMS 75  | 1220 x 610 | 37 | 21737     |
| Laminar Flow Module FMS 93  | 1525 x 610 | 45 | 21739     |
| Laminar Flow Module FMS 112 | 1830 x 610 | 50 | 21741     |

(Height of module: 420 mm)

#### Laminar Flow Box FBS-Series – Complete Laminar Flow Box incl. filter module

##### Dimensions Laminar Flow Box - FBS-Series:

| Type                     | (W x D mm) | Order No. |
|--------------------------|------------|-----------|
| Laminar Flow Box FBS 24  | 610 x 400  | 21732     |
| Laminar Flow Box FBS 37  | 610 x 610  | 21734     |
| Laminar Flow Box FBS 56  | 915 x 610  | 21736     |
| Laminar Flow Box FBS 75  | 1220 x 610 | 21738     |
| Laminar Flow Box FBS 93  | 1525 x 610 | 21740     |
| Laminar Flow Box FBS 112 | 1830 x 610 | 21742     |

(Height inside: 700 mm, total height: 1200 mm)

### Equipment Versions:

- Version 1:** stainless steel base frame on castors - for all versions
- Version 2:** customised laboratory cabinet mounted upon request
- Version 3:** Polycarbonate version, UV-resistant
- Version 4:** interior lighting

- Version 5:** interior lighting with UV-light
- Version 6:** 230 V connection inside, triple wall socket
- Version 7:** automatic re-adjustment of air flow with filter exchange indicator
- Version 8:** stainless steel ceiling mount