

Accessories - Your daily assistance, practical and safe



PE / PP trays

- Increased protection of drawers when storing mildly corrosive liquids and in case of any spillages
- Removable; safe disposal of leakage and easy cleaning

Rubber mats

- For optimal stability of containers on the drawer




PE trays for shelves

- Increased protection of shelves when storing mildly corrosive liquids and in case of any leakages
- Removable; safe disposal of leakage and easy cleaning

Bottle racks

- Help to store small bottles clearly arranged on the drawer



Box for chemicals

- Separate storage of mildly corrosive chemicals inside the cabinet
- Optimal for safe indoor transport

Storage box

- Lockable, separate storage of hazardous materials inside the cabinet



Fire protected add-on

- Allows the installation of pipes or cables through the side or rear walls of Typ 90 cabinets, without compromising the fire resistance according to EN 14470-1
- Certificate by an independent material testing institute available upon request