

# «« Safety Storage Cabinets



## Type 30 Underbench cabinets



Model F30.60.140-S with 2 drawers



Model F30.60.110-S with 1 drawer



Extraction air connection $\varnothing$	mm	75
Change of air (10 times)	m <sup>3</sup> /h	3
Pressure drop	Pa	1

Extraction air connection $\varnothing$	mm	75
Change of air (10 times)	m <sup>3</sup> /h	2
Pressure drop	Pa	1

### Safety Storage Cabinets

**30 minutes** fire resistance in accordance with **EN 14470-1 (30 minutes fire rated - Type 30)**

GS and CE conformity

Labelling in accordance with **EN 14470-1**

### F30.60.140-S

Dimensions W x D x H (mm):

**exterior approx. righth: 1400 x 550 x 600**

**interior approx. left: 770 x 445 x 510**

**approx. 470 x 445 x 510**

Weight: **approx. 145 kg**

### Drawers

Dimensions W x D x H (mm):

Left: **approx. 722 x 442 x 80**

Capacity: **approx. 22 ltr.**

Right: **approx. 422 x 442 x 80**

Capacity: **approx. 12 ltr.**

Load capacity (evenly distributed load): **2 x 50 kg**

Useable interior height:

**approx. 470 mm** (without perforated metal plate insert)

### F30.60.110-S

Dimensions W x D x H (mm):

**exterior approx. 1100 x 550 x 600**

**interior approx. 980 x 445 x 510**

Weight: **approx. 120 kg**

### Drawer

Dimensions W x D x H (mm):

**approx. 930 x 442 x 80**

Capacity: **approx. 27 ltr.**

Load capacity (evenly distributed load): **50 kg**

Useable interior height:

**approx. 470 mm** (without perforated metal plate insert)

### Type 30 Underbench cabinets with 1 or 2 drawers

Model	Material, colour	recommended quantity	Order Code
<b>F30.60.140-S</b> , Underbench cabinet, with drawers	light grey, RAL 7035		<b>20241</b>
<b>Interior equipment</b>			
Second level drawer, left side	sheet-steel, plastic laminated		<b>9609</b>
Second level drawer, right side	sheet-steel, plastic laminated		<b>9608</b>
Movable plinth, height 30 mm	sheet-steel, plastic laminated		<b>21714</b>
Movable plinth, height 65 mm (30 mm + 35 mm faceplate)	sheet-steel, plastic laminated		<b>21896</b>
<b>F30.60.110-S</b> , Underbench cabinet, with drawer	light grey, RAL 7035		<b>20237</b>
<b>Interior equipment</b>			
Second level drawer	sheet-steel, plastic laminated		<b>8053</b>
Movable plinth, height 30 mm	sheet-steel, plastic laminated		<b>19705</b>
Movable plinth, height 65 mm (30 mm + 35 mm faceplate)	sheet-steel, plastic laminated		<b>21895</b>


## Type 30 Underbench cabinets



Model F30.60.110-2S with 2 drawers



Model F30.60.110-T with 2 wing doors



Extraction air connection $\varnothing$	mm	75
Change of air (10 times)	m <sup>3</sup> /h	2
Pressure drop	Pa	1

### Safety Storage Cabinets

**30 minutes** fire resistance in accordance with **EN 14470-1 (30 minutes fire rated - Type 30)**

GS and CE conformity

Labelling in accordance with **EN 14470-1**

### F30.60.110-2S

Dimensions W x D x H (mm):

**exterior approx. 1100 x 550 x 600**

**interior 2 x approx. 470 x 445 x 510**

Weight: **approx. 120 kg**

### Drawers

Dimensions W x D x H (mm):

**2 x approx. 422 x 442 x 80**

Capacity: **approx. 22 ltr.**

Load capacity (evenly distributed load): **2 x 50 kg**

Useable interior height:

**approx. 470 mm** (without perforated metal plate insert)

### F30.60.110-T

Dimensions W x D x H (mm):

**exterior approx. 1100 x 550 x 600**

**interior approx. 980 x 445 x 500**

Weight: **approx. 110 kg**

### Shelf

Dimensions W x D x H (mm):

**approx. 970 x 415 x 20**

Load capacity (evenly distributed load): **75 kg**

Capacity of the bottom collecting sump: **approx. 30 ltr.**

Useable interior height:

**approx. 460 mm** (without perforated metal plate insert)

### Type 30 Underbench cabinets with 2 drawers / 2 wing doors

Model	Material, colour	recommended quantity	Order Code
F30.60.110-2S, Underbench cabinet, with drawers	light grey, RAL 7035		<b>20235</b>
<b>Interior equipment</b>			
Second level drawer, left or right	sheet-steel, plastic laminated		<b>9608</b>
F30.60.110-T, Underbench cabinet, with wing doors	light grey, RAL 7035		<b>20239</b>
<b>Interior equipment</b>			
Bottom collecting sump (V=20 ltr.)	sheet-steel, plastic laminated	1	<b>8372</b>
Perforated insert	sheet-steel, plastic laminated	1	<b>8374</b>
Shelf	sheet-steel, plastic laminated	1	<b>8376</b>
<b>Accessories (both models)</b>			
Movable plinth, height 30 mm	sheet-steel, plastic laminated		<b>19705</b>
Movable plinth, height 65 mm (30 mm + 35 mm faceplate)	sheet-steel, plastic laminated		<b>21895</b>