

Type 90 Safety Underbench Cabinets with new SOFT-CLOSE Technology

Gründau, 1st March 2016 – Effortless and quiet: the asecos soft-closing technology ensures that the drawers of the safety underbench cabinets are securely closed. The SOFT-CLOSE system additionally helps to avoid potential damage of stored hazardous materials when the doors are pushed too hard.

asecos Type 90 safety underbench cabinets are GS certified and offer a maximum fire resistance of 90 minutes. Available in 14 variants with four widths and two depths 60cm and now 80cm.

The exclusive feature of our new safety storage underbench cabinet is the special height of 80cm offering storage of containers up to 33cm high on two levels. It is possible to store up to 32 Winchesters dependent on size. To prevent unauthorized use doors and drawers are lockable.

For further information, please visit: <http://www.asecos.co.uk>

Contact person:

asecos GmbH
Christina Bien
Tel. +49 6051 9220-523
Fax: +49 6051 9220-727
E-Mail: c.bien@asecos.com

asecos GmbH

asecos is the leading German manufacturer of highly-insulated safety storage cabinets for flammables. For more than 20 years asecos is actively involved in international associations and working groups to develop innovative and sustainable products to meet the challenges of storing hazardous materials in 21st century laboratories and industrial facilities. asecos cabinets are built according to European and Ameri-

PRESS RELEASE



2/2

can standards and can be integrated into casework or used as free-standing storage solutions. In addition asecos also offers many products for the safe handling of chemicals as well as a variety of extraction systems. Abroad asecos is represented by own subsidiaries in the Netherlands, France, Spain and the UK and maintains a global network of specialized dealers and distributors.

For further information, please visit: www.asecos.com