

Everything is clear! with the asecos Q-DISPLAY-30

Gründau, 25th of July 2016 – The Q-DISPLAY-30 fire resistant safety cabinet with an extra-large glass insert allows you to see at a glance all materials stored within. This safety storage cabinet is ideal for installation in laboratories, display or dispensing points and retail.

The extremely robust fire resistant Q-DISPLAY-30 cabinet is manufactured to European Standard EN 14470-1. At 86cm wide there is room for plenty of Chemicals. The triple hinged door with scratch and impact resistant surfaces, remains open in any position and is self-closing in the event of fire. Furthermore, the cabinet is equipped with a secure cylinder lock and an advanced ventilation system ready for connection to external LEV or asecos recirculating air filtration systems.

The integrated transport base allows the cabinet to be easily relocated, in the event of the re-tasking of laboratories, deep cleaning and company relocations. In addition, a wide range of interior equipment is available: Drawers, height-adjustable shelves, bottom collecting sump for safe collection of spillages and even metal-free interior equipment for the storage of aggressive hazardous materials.

For further information, please visit: <http://www.asecos.co.uk>

Contact person:

asecos GmbH

Christina Bien

Tel. +49 6051 9220-523

Fax: +49 6051 9220-727

E-Mail: c.bien@asecos.com

PRESS RELEASE



2/2

asecos GmbH

asecos is the leading German manufacturer of highly-insulated safety storage cabinets for flammables. For more than 20 years asecos is actively involved in international associations and working groups to develop innovative and sustainable products to meet the challenges of storing hazardous materials in 21st century laboratories and industrial facilities. asecos cabinets are built according to European and American standards and can be integrated into casework or used as free-standing storage solutions. In addition asecos also offers many products for the safe handling of chemicals as well as a variety of extraction systems. Abroad asecos is represented by own subsidiaries in the Netherlands, France, Spain and the UK and maintains a global network of specialized dealers and distributors.

For further information, please visit: www.asecos.com